



18962

## Phenyl Isothiocyanate (PITC) extrapure AR, 99%

Part B

CAS : 103-72-0

Molecular Formula : C<sub>7</sub>H<sub>5</sub>NS

Molecular Weight : 135.19

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Almost colourless to yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.130-1.132
Refractive Index (20°C)	1.649-1.650
Water (KF)	max. 0.1%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2927 , IMCO Class No.:6.1(8) , Packing Group:II
HSN Code	
25 ml	29309099 (GST 18%)
5 x5 ml	29309099 (GST 18%)
100 ml	29309099 (GST 18%)

### Available Packages

25 ml

100 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com