



20139

Trimethyl Chlorosilane (TMCS) extrapure, 98%

Part B

CAS : 75-77-4

Molecular Formula : C₃H₉ClSi

Molecular Weight : 108.64

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay	min. 98%
Density (g/ml) @ 20°C	0.858-0.860
Refractive Index (20°C)	1.388-1.389
Boiling Range	56-58°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1298 , IMCO Class No.:3(8) , Packing Group:II
HSN Code	
1000 ml	29319090 (GST 18%)
25 ml	29319090 (GST 18%)
250 ml	29319090 (GST 18%)
2500 ml	29319090 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle
25 ml	Glass Bottle
250 ml	Glass Bottle
2500 ml	Glass Bottle

Available Packages

250 ml

1000 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com