



24564

## Polyvinylpyrrolidone K90 (Povidone, PVP K-90)

Part A

CAS : 9003-39-8

Molecular Formula : (C<sub>6</sub>H<sub>9</sub>NO)<sub>n</sub>

Molecular Weight : ~360000

### General Information

Storage	25 to 40°C (Room Temperature)
HSN Code	
100 Gms	39059990 (GST 18%)
500 Gms	39059990 (GST 18%)
5 Kg	39059990 (GST 18%)

### Available Packages

100 Gms  
500 Gms  
5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com