



28065

Styrene Monomer (Stab w/ tBC) extapure, 99.5%

Part B

CAS : 100-42-5

Molecular Formula : C₆H₅CH=CH₂

Molecular Weight : 104.15

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to yellow oily
Appearance (Form)	Liquid
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.906-0.908
Refractive Index (20°C)	1.544-1.546
Boiling Point	145-146°C
Stabilizer (4-tert-butyl catechol)	10-15ppm

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2056 , IMCO Class No.:3 , Packing Group:II
HSN Code	
500 ml	29025000 (GST 18%)
2500 ml	29025000 (GST 18%)

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com