



32385

Carrageenan (Kappa) Refined ex. Eucheuma Cottonii

Part A

CAS : 11114-20-8

Molecular Formula : C₂₄H₃₆O₂₅S₂-2

Specifications

Appearance (Colour)	White to off-white to pale yellow
Appearance (Form)	Powder
Loss on drying	max.12%
Ash	max. 40%
Potassium gel strength (1.5% in 0.3% aq.KCl soln.)	min. 1300 g/cm ²
Viscosity (0.3% aq. at 25°C)	10 - 40 mPas

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	13023900 (GST 18%)
25 Gms	13023900 (GST 18%)
500 Gms	13023900 (GST 18%)
Type of Packing	
25 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com