



36805

## Benzoic Acid pure, 99%

Part A

CAS : 65-85-0

Molecular Formula : C<sub>7</sub>H<sub>6</sub>O<sub>2</sub>

Molecular Weight : 122.12

### Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Compound or crystal or flakes
Solubility (Turbidity) 10% solution in methanol	Clear
Assay (by Neutralization)	min. 99.5%
Melting Point	121 - 122°C
Sulphated Ash	max. 0.05%
Total chlorine	max. 0.1%
Heavy Metals (Pb)	max. 0.0005%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	29163110 (GST 18%)
500 Gms	29163110 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms

5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com