



37020

Azomethine H Monosodium Salt Hydrate, 97%

Part B

CAS : 206752-32-1

Molecular Formula : C₁₇H₁₂NO₈S₂Na.XH₂O

Molecular Weight : 445.40

Specifications

Appearance (Colour)	Yellow to orange
Appearance (Form)	Powder
Assay (On dry basis)	min. 97%
Absorbance (A) of 1% solution (pH 5.1) in a 1cm cell @236nm	min. 1200

Other Information

Application

Azomethine H is a colorimetric reagent for boron detection.

It forms an orange complex with boron in aqueous solution. The molar absorptivity of this complex is 1.01x10⁴ at 415nm. The detection range of boron in sample solutions is 1-6µm per ml or 1.0-6ppm. To detect boron in plant samples EDTA is used to mask copper, iron and aluminium ions.

Boron determinations are used in,
1. Soil, plants, composts, manure, water and nutrient solutions
2. Glass and Steel samples (microgram levels of boron)

Azomethine H is also used in Electrocyclization Reactions in the synthesis of Martinelllic acid, Spirotryprostatin A, Benzodiazepinones.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29222190 (GST 18%)
25 Gms	29222190 (GST 18%)
5 Gms	29222190 (GST 18%)

Type of Packing

1 Gms

Glass Bottle

25 Gms

Glass Bottle

5 Gms

Glass Bottle

Available Packages

1 Gms

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
