



38063

Glycolic Acid extrapure, 70% soln in water

Part A

CAS : 79-14-1

Molecular Formula : C₂H₄O₃

Molecular Weight : 76.05

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Pale yellow
Appearance (Form)	Liquid
Assay (NaOH Titration)	min. 70%
Density (g/ml) @ 20°C	~1.25 g/ml

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3261 , IMCO Class No.:8 , Packing Group:II
HSN Code	
500 ml	29181190 (GST 18%)
2500 ml	29181190 (GST 18%)

Available Packages

500 ml
2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com