



41030

## Polyvinylpyrrolidone K30 (Povidone, PVP K-30)

Part A

CAS : 9003-39-8

Molecular Formula : (C<sub>6</sub>H<sub>9</sub>NO)<sub>x</sub>

Molecular Weight : ~40000

### Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Nitrogen (N) content	11.5 - 12.5%
pH (5% aq. solution)	3.0 - 7.0
K-value	~ 30

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	39059990 (GST 18%)
25 Kg	39059990 (GST 18%)
5 Kg	39059990 (GST 18%)
500 Gms	39059990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

100 Gms  
500 Gms  
5 Kg  
25 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com