



46651

6-Benzyladenine (6-BAP) extrapure AR, 99%

Part B

CAS : 1214-39-7

Molecular Formula : C₁₂H₁₁N₅

Molecular Weight : 225.25

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 1% solution in 0.1N HCl	Clear
Solubility (Colour) 1% solution in 0.1N HCl	Colourless
Assay (NT)	min. 99%
Loss on drying	max. 0.1%
Sulphated Ash	max. 0.05%
Melting Point	230 - 233°C

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	98020000 (GST 18%)
25 Gms	98020000 (GST 18%)
5 Gms	98020000 (GST 18%)
Type of Packing	
5 Gms	Glass Bottle

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com