



48088

Charcoal Activated (950) Granular extrapure

Part A

CAS : 7440-44-0

Molecular Weight : 12.01

Specifications

Appearance (Colour)	Black
Appearance (Form)	Granular compound
pH (2 % aq. suspension)	6.0 - 9.0
Iodine value	min. 950 mg/gm
Loss on drying	max 10%
Sulphated Ash	max 4%
Chloride (Cl)	max 0.1%
Sulphate (SO ₄)	max 0.05%
Iron (Fe)	max 0.05%

Other Information

Notes	(Forbidden by Air)
-------	--------------------

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1362 , IMCO Class No.:4.2 , Packing Group:III
HSN Code	
5 Kg	38029019 (GST 18%)
500 Gms	38029019 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

500 Gms
5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com