



51562

Chlorogenic Acid Hemihydrate extrapure, 98%

Part B

CAS : 6001-76-9

Molecular Formula : C₁₆H₁₈O₉.1/2H₂O

Molecular Weight : 363.32

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% solution in Methanol	Clear
Assay (NaOH titration)	min. 98.0%
Melting Point	205 - 215°C

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	98020000 (GST 18%)
100 Mg	98020000 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle

Available Packages

100 Mg

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com