



52392

## L-Arabinose extrapure

Part B

CAS : 87-72-9, 5328-37-0

Molecular Formula : C<sub>5</sub>H<sub>10</sub>O<sub>5</sub>

Molecular Weight : 150.13

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Melting Point	154° - 160° C
Specific rotation [?]	+103° to +106° (C = 5, water)
Sulphated Ash	max. 0.1%
Water (KF)	max. 0.5%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29400000 (GST 18%)
25 Gms	29400000 (GST 18%)
500 Gms	29400000 (GST 18%)
Type of Packing	
10 Gms	Plastic Bottle
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

### Available Packages

100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com