



55664

Tris(2-carboxyethyl) Phosphine Hydrochloride (TCEP) extrapure AR, 98%

Part B

CAS : 51805-45-9

Molecular Formula : C₉H₁₅O₆P.HCl

Molecular Weight : 286.65

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay (Iodometric)	min. 98%

Other Information

Notes	A superior reducing agent than DTT.
-------	-------------------------------------

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:3261 , IMCO Class No.:8 , Packing Group:II
HSN Code	
1 Gms	29211190 (GST 18%)
25 Gms	29211190 (GST 18%)
5 Gms	29211190 (GST 18%)
Type of Packing	
5 Gms	Glass Bottle

Available Packages

1 Gms
5 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com