



55895

6x TriColor Loading Dye

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Dark Blue
Electrophoresis test	Dye in 1X Concentration is mixed with Plasmid DNA and the mobility is tested by agarose gel electrophoresis
Stability	It should be stable when subjected to repeated freeze-thaw

Other Information

Description

6X Tri-colour Loading Dye is a pre-mixed loading buffer with three tracking dyes for agarose and non-denaturing polyacrylamide gel electrophoresis. EDTA is included to chelate magnesium (up to 10 mM) in enzymatic reactions, thereby stopping the reaction. Bromophenol Blue, xylene cyanol and orange G are the standard tracking dyes for electrophoresis. Analysis of large DNA molecules. Preparation of DNA ladders, markers and samples for loading on agarose or polyacrylamide gels

Application

Features:
Tri-color tracking of DNA migration during electrophoresis. No DNA masking during gel exposure to UV light. EDTA binds divalent metal ions and inhibits metal dependent nucleases.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 x1 ml	38229090 (GST 12%)

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com