



58389

## Heparin Sodium Salt ex. Bovine Intestinal Mucosa, 150 IU/mg

Part B

CAS : 9041-08-1

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% aq. solution	Clear
Activity	min. 150 IU/mg
pH (1% aq. solution)	5.0 -7.5
Water (KF)	max. 8%
Sulphated Ash	28 - 41%
Heavy Metals (Pb)	max. 0.003%
Absorbance (A) of 0.4% aq. solution in a 1cm cell v/s H <sub>2</sub> O	
@260nm	max. 0.20
@280nm	max. 0.15

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 K.Units	38229090 (GST 12%)
20 K.Units	38229090 (GST 12%)
Type of Packing	
100000 Units	Glass Bottle
2.5 M.Units	Glass Bottle
20000 Units	Glass Bottle

### Available Packages

20 K.Units

100 K.Units

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com