



60020

Z-L-Arginine extrapure, 99%

Part B

CAS : 1234-35-1

Molecular Formula : C₁₄H₂₀N₄O₄

Molecular Weight : 308.3

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Assay (NT) (On dry basis)	min. 99%
Melting Point	171° - 174°C
Water (KF)	max. 4%

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
HSN Code	
25 Gms	29379020 (GST 18%)

Available Packages

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com