



64236

Xanthan Gum ExiPlus, Multi-Compendial

Part A

CAS : 11138-66-2

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	White to yellow to beige
Appearance (Form)	Free flowing powder
pH (1% aq. solution)	6.0 - 8.0
Loss on drying	max. 15%
Ash	6.5-16%
Heavy Metals (Pb)	max. 0.003%
Pyruvic acid	min.1.5%
Viscosity (1% in 1% KCL solution. at 25°C)	min. 600 mpas
Particle Size	Min 92% from through Mesh

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	39139090 (GST 18%)
500 Gms	39139090 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

Available Packages

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com