



64407

## TMB Solution A (1% TMB in DMSO)

Part B

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellowish orange
Appearance (Form)	Solution
Composition	1% 3,3,5,5-TMB in DMSO
Suitability as chromogenic substrate	Passes

### Other Information

Application

Application:  
To prepare a working solution for ELISA, dilute 1.0 mL of TMB solution A with 9.0 mL of phosphate-citrate buffer pH 5.2, and add 2.0 µL of TMB solution B. This solution is now ready to be used for ELISA involving an HRP-conjugated Ab.

Suitable for ELISA and HRP (Peroxidase) related applications

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 ml	38229090 (GST 12%)
Type of Packing	
100 ml	Glass Bottle

### Available Packages

100 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com