



65240

NTP Set Ultrapure (100mM each)

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Colourless
Appearance (Clarity)	Clear
pH	8 ± 0.2
Absence of nucleases	None detected

Other Information

Description

The set consists of 100mM aqueous solutions of ATP, UTP, CTP and GTP at pH 7.5 as sodium salts each supplied in a separate vial. Since the nucleotides are provided separately, the NTP Set offers maximum flexibility in preparation of reaction mixes for different applications. The set is free of Dnase, Rnase, Protease and Nicking activity.

Applications

In- vitro transcription reactions.
Production of RNA probes.

Includes

4 x 0.25ml
ATP Solution Ultrapure (100mM) –
0.25ml
CTP Solution Ultrapure (100mM) –
0.25ml
GTP Solution Ultrapure (100mM) –
0.25ml
UTP Solution Ultrapure (100mM) –
0.25ml

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
4 x0.25 ml	38229090 (GST 12%)

Available Packages

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com