



67577

## Antimony Metal Ingots ultrapure, 99.999%

Part A

CAS : 7440-36-0

Molecular Formula : Sb

Molecular Weight : 121.75

### Specifications

Appearance (Colour)	Shiny metal
Appearance (Form)	Ingots
Assay	min. 99.999%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2871 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
10 Gms	81109000 (GST 18%)
25 Gms	81109000 (GST 18%)
Type of Packing	
25 Gms	Plastic Bottle

### Available Packages

10 Gms  
25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com