



67883

Trimethyl Chlorosilane (TMCS) 1M Solution in THF

Part B

CAS : 75-77-4

Molecular Formula : C₃H₉ClSi

Molecular Weight : 108.64

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Molarity (Argentometric titration)	0.95-1.1M

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	24 Months
HSN Code	
100 ml	29319090 (GST 18%)

Available Packages

100 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com