



68913

meta-Topolin (mT) for molecular biology, 98%

Part B

CAS : 75737-38-1

Molecular Formula : C₁₂H₁₁N₅O

Molecular Weight : 241.25

Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Assay	min. 98%
DNase, RNase, Protease	Not detected

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Mg	98020000 (GST 18%)
25 Mg	98020000 (GST 18%)
Type of Packing	
100 Mg	Glass Bottle

Available Packages

25 Mg

100 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com