



69114

p-Nitroaniline (High Purity) extrapure AR, 99.5%

Part A

CAS : 100-01-6

Molecular Formula : C₆H₆N₂O₂

Molecular Weight : 138.13

Specifications

Appearance (Colour)	Yellow (may turn dark yellow to brown on storage)
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 3% solution in Methanol	Clear
Solubility (Colour) 3% solution in Methanol	Yellow
Assay (GC)	min. 99.5%
Melting Point	147 - 149°C
Sulphated Ash	max. 0.05%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%
Water (KF)	max. 0.5%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1661 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
100 Gms	29214226 (GST 18%)
25 Gms	29214226 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com