



70984

## CFC Supplement

Part D

### Specifications

Appearance (Colour)	White, fine powder
Solubility	Soluble in distilled water
Cultural Response	Inoculate and incubate at $25 \pm 2^\circ\text{C}$ for 24-48 hours. (Aseptically add the rehydrated contents of 1 vial as per direction) Growth observed on Pseudomonas Agar Base (I) (90912).
Organism	Pseudomonas aeruginosa ATCC 27853
Inoculum (cfu) 10-100	Growth : Inhibited
Organism	Staphylococcus aureus ATCC 25923
Inoculum	Growth : Inhibited
Organism	Proteus vulgaris ATCC 13315
Inoculum	Growth : Inhibited

### Other Information

#### Applications

Used for selective isolation of Pseudomonas species.

#### Composition

Cetrimide	5.00 mg/vial
Fucidin	5.00 mg/vial
Cephaloridine	25.00 mg/vial

#### Directions

1. Aseptically add 2 ml of sterile distilled water to 1 vial.
2. Mix well and aseptically add to sterile cooled media at  $45-50^\circ\text{C}$  and pour in petriplates.
3. Add 2ml (1 vial) of supplement to 500 ml of Pseudomonas Agar Base (I) (90912).

### General Information

Storage	2 to $8^\circ\text{C}$ (Refrigerate)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 vl	38229090 (GST 12%)
5 x5 vl	38229090 (GST 12%)

## Available Packages

5 vI

5 x5 vI

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---