



72866

Molecular Sieve 4A (0.4nm) (1.6-3mm Beads)

Part A

CAS : 70955-01-0

Molecular Formula : $\text{Na}_{12}[(\text{AlO}_2)_{12}(\text{SiO}_2)_{12}]\cdot x\text{H}_2\text{O}$

Specifications

Appearance (Colour)	Pale yellow to brown
Appearance (Form)	Beads
Bed crushing strength	min. 85%
Water absorption capacity	min. 20%
Bulk density	~740 kg/m ³
Pore diameter	4 Å or 0.4 nm
Loss on ignition	max. 3%
Particle Size	6 - 12 mesh
Particle diameter	1.6 - 3 mm

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	28399090 (GST 18%)
250 Gms	28399090 (GST 18%)
5 Kg	28399090 (GST 18%)

Available Packages

250 Gms
1 Kg
5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com