



73317

## Epichlorohydrin extrapure AR, 99%

Part A

CAS : 106-89-8

Molecular Formula : C<sub>3</sub>H<sub>5</sub>ClO

Molecular Weight : 92.53

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.180-1.183
Refractive Index (20°C)	1.437-1.438
Boiling Range	115-117°C
Non Volatile Matter	max. 0.005%
Acidity (CH <sub>3</sub> COOH)	max. 0.001%
Water (KF)	max. 0.1%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2023 , IMCO Class No.:6.1(3) , Packing Group:II
HSN Code	
2500 ml	29103000 (GST 18%)
500 ml	29103000 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle

### Available Packages

500 ml

2500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com