



83826

## DEAE-Dextran 500 ex. Leuconostoc Sp.

Part B

CAS : 9015-73-0

### Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Powder
Average Molecular weight	made from Dextran 400 -750 KDa
Nitrogen (N) content	2.5 - 4.5%
Loss on drying	max. 10%

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
HSN Code	
1 Gms	39140090 (GST 18%)
5 Gms	39140090 (GST 18%)

### Available Packages

1 Gms  
5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com