



85441

## D-Amygdalin extrapure, 97%

Part B

CAS : 29883-15-6

Molecular Formula : C<sub>20</sub>H<sub>27</sub>NO<sub>11</sub>

Molecular Weight : 457.44

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Assay	min. 97%
Specific rotation [?]	-35° to -45° (C=2, Water)
Melting Point	210° - 215°C
Loss on drying	max. 5%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.0005%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29269000 (GST 18%)
10 Gms	29269000 (GST 18%)
25 Gms	29269000 (GST 18%)
Type of Packing	
10 Gms	Plastic Bottle
25 Gms	Plastic Bottle

### Available Packages

1 Gms

10 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com