



86320

## 3-Cyanophenylboronic Acid extrapure, 97%

Part B

CAS : 150255-96-2

Molecular Formula : C<sub>7</sub>H<sub>6</sub>BNO<sub>2</sub>

Molecular Weight : 146.9

### Specifications

Appearance (Colour)	Off - white
Appearance (Form)	Crystalline powder
Assay (NaOH titration)	min. 97.0%

### General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29319090 (GST 18%)
25 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle
5 Gms	Glass Bottle

### Available Packages

1 Gms

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com