



88109

Isomaltotetraose extrapure, 95%

Part B

CAS : 35997-20-7

Molecular Formula : C₂₄H₄₂O₂₁

Molecular Weight : 666.58

Other Information

Application

Isomaltotetraose, can be used for immunization of hybridoma antibodies. the use of synthetic glycolipids as antigens enables one to study the gene control of antibody responses to glycolipids and to investigate the combining-site specificities of antibodies to a single antigenic determinant.

General Information

Storage

8 to 25°C (Cool & Dry Area)

HSN Code

5 Mg

29400000 (GST 18%)

Available Packages

5 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com