



89465

Florisil 60-100 mesh for chromatography

Part B

CAS : 1343-88-0

Molecular Formula : MgO3Si

Molecular Weight : 100.39

Specifications

| | |
|----------------------|--|
| Appearance (Colour) | White |
| Appearance (Form) | Free flowing Fine powder, may contain some dark particle |
| Chloride (Cl) | max. 0.02% |
| Iron (Fe) | max. 0.02% |
| Mesh size (60 - 100) | min. 80% |

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 1 Kg | 38249900 (GST 18%) |
| 100 Gms | 38249900 (GST 18%) |
| 250 Gms | 38249900 (GST 18%) |
| Type of Packing | |
| 250 Gms | Plastic Bottle |

Available Packages

100 Gms

250 Gms

1 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com