



90739

Uric Acid extrapure AR, 99%

Part B

CAS : 69-93-2

Molecular Formula : C₅H₄N₄O₃

Molecular Weight : 168.11

Specifications

Appearance (Colour)	White to off - white
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 5% solution in hot 1N NaOH	Clear
Solubility (Colour) 5% solution in hot 1N NaOH	Colourless to pale yellow
Assay (T)	min. 99%
Sulphated Ash	max. 0.1%
Water (KF)	max. 0.5%

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29335990 (GST 18%)
25 Gms	29335990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com