



92300

## Benzyl Alcohol extrapure, 99%

Part A

CAS : 100-51-6

Molecular Formula : C<sub>7</sub>H<sub>8</sub>O

Molecular Weight : 108.14

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.044 - 1.048
Refractive Index (20°C)	1.539 - 1.541
Boiling Range	202 - 205°C
Non Volatile Matter	max. 0.05%
Acidity (C <sub>6</sub> H <sub>5</sub> COOH)	max. 0.025%
Chlorinated compounds (Cl)	max. 0.02%
Miscibility with water (1:30)	Clear & colourless
Water (KF)	max. 0.2%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1000 ml	29062100 (GST 18%)
25 Ltr	29062100 (GST 18%)
2500 ml	29062100 (GST 18%)
500 ml	29062100 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

## Available Packages

500 ml

1000 ml

2500 ml

25 Ltr

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---