



92604

## Cyclohexane Dried, 99.5%, water 0.005%

Part A

CAS : 110-82-7

Molecular Formula : C<sub>6</sub>H<sub>12</sub>

Molecular Weight : 84.16

### Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.777-0.779
Refractive Index (20°C)	1.425-1.426
Boiling Range	80-81°C
Freezing Point	6-7°C
Non Volatile Matter	max. 0.001%
Acidity (HCl)	max. 0.001%
Substance Discoloured by H <sub>2</sub> SO <sub>4</sub>	Passes test
Water (KF)	max. 0.005%

### Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1145 , IMCO Class No.:3 , Packing Group:II
HSN Code	
1000 ml	29021100 (GST 18%)
500 ml	29021100 (GST 18%)
Type of Packing	
1000 ml	Glass Bottle

## Available Packages

500 ml

1000 ml

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---