



93633

## UTP Solution Ultrapure (100mM)

Part E

### Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Colourless
Appearance (Clarity)	Clear
Absence of nucleases	DNase, RNase, Protease not detected
pH	8.2± 0.2

### Other Information

Description

This is a vial of 100mM aqueous solution of UTP (uridine 5'-triphosphate), at pH 7.5 as sodium salt. The preparation is free of Dnase, Rnase, Protease and Nicking activity

#### Applications

In vitro transcription  
aRNA synthesis  
siRNA synthesis  
RNA amplification

Includes

UTP Solution Ultrapure (100mM) –  
0.25ml

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
0.25 ml	38229090 (GST 12%)

### Available Packages

0.25 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com