



98669

n-Hexanol extrapure AR, 99%

Part B

CAS : 111-27-3

Molecular Formula : C₆H₁₄O

Molecular Weight : 102.18

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay (GC)	min. 99%
Single impurity	max. 0.4%
Density (g/ml) @ 20°C	0.817-0.819
Refractive Index (20°C)	1.418-1.419
Boiling Range	156-157°C
Non Volatile Matter	max. 0.001%
Acidity (CH ₃ COOH)	max. 0.001%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2282 , IMCO Class No.:3 , Packing Group:III
HSN Code	
250 ml	29051690 (GST 18%)
2500 ml	29051690 (GST 18%)
500 ml	29051690 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

250 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com