



98680

## Bathophenanthroline Disulphonate Disodium Salt Hydrate extrapure, 99%

Part B

CAS : 52746-49-3, 98645-86-4, 53744-

Molecular Formula : C<sub>24</sub>H<sub>14</sub>N<sub>2</sub>Na<sub>2</sub>O<sub>6</sub>S<sub>2</sub>.xH<sub>2</sub>O

Molecular Weight : 536.49 (anhydrous basis)

### Specifications

Appearance (Colour)	White to Orange and Faint Beige to Beige and Pink
Appearance (Form)	Powder
Assay (HPLC) (on anhydrous basis)	min. 99%
Identification (IR Spectrum)	Passes test

### Other Information

Application	A water soluble reagent for determination of Iron (Fe).
-------------	---

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29335990 (GST 18%)
100 Mg	29335990 (GST 18%)
250 Mg	29335990 (GST 18%)

### Available Packages

250 Mg

1 Gms

5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com